


## IS THE HAIRMAX® RIGHT FOR YOU?

HairMax® is most effective in women and men in early to moderate stages of hair loss and thinning. If your hair loss looks like the images below, HairMax is recommended for you.

### FEMALE



### MALE



## DOCTOR RECOMMENDED

Physicians worldwide recommend HairMax to their patients with hereditary hair loss. HairMax is supported by a Medical Advisory Board of established key authoritative figures in the fields of dermatology and hair restoration.

## FREQUENTLY ASKED QUESTIONS

### How often do I need to use my LaserBand 82?

It is recommended to use your LaserBand three times a week for as little as 90 seconds, with at least a day between each treatment.

### Will results improve if I use my LaserBand 82 more than 3 times a week?

Three times a week has been found to be the optimal treatment time. Using your LaserBand more than 3 times a week has not shown any additional benefits.

### How long before I see results?

Most users start to see benefits in 12 - 16 weeks.

### How long will I need to continue hair loss treatments with the LaserBand 82?

To keep seeing benefits use your LaserBand as directed on an ongoing basis. If you do not use the device as recommended, there is a strong possibility that you will not see benefits or your hair loss will regress to where you were prior to using the LaserBand.

### Can I use my LaserBand 82 on wet or dry hair?

You can use your device on dry or damp hair. It is important that your scalp be clean and free of styling products. A build-up of topicals such as oils, gels, sprays, etc. can block the laser light from effectively reaching the hair follicles.

### Do I have to use any special shampoos or products with my LaserBand 82?

No, however, we recommend using the HairMax for den•si•ty collection. These sulfate-free, low pH products are gentle on your thinning hair and complement the laser treatment.



Lexington International, LLC  
1040 Holland Drive, Boca Raton, FL 33487  
Tel. +1.561.417.0200 | U.S. Toll Free 1.800.9REGROW  
FAX. +1.561.892.0747

[www.HairMax.com](http://www.HairMax.com)

LIT-CONBROLEX-LB82 03 07/17 © Copyright 2017 Lexington Intl., LLC. All rights reserved.



## BREAKTHROUGH TREATMENT FOR HAIR LOSS & THINNING HAIR

HairMax® LaserBand 82




DESIGN BY *pininfarina*

**HAIRMAX®**  
The science of hair growth.™

## HOW LASER LIGHT THERAPY PROMOTES HAIR GROWTH

HairMax laser treatments work through the process of photobiostimulation, delivering low level laser light to your hair follicles to activate them for stronger hair growth.




This boost of nourishing light energizes hair follicles at a cellular level to stimulate growth factors, reverse hair thinning, restore your hair's natural growth cycle and help your hair to grow fuller, denser, longer and stronger!

## FASTEST TREATMENT AVAILABLE





# LaserBand 82 - The Ultimate Hair Growth Laser Device



82 medical grade lasers over 3 areas of the scalp provide full coverage equating to 246 lasers

NO LEDs

DESIGN BY *pininfarina*

- Stimulates Hair Growth
- Treats Hair Loss
- Reverses Thinning
- Increases Density & Fullness
- Revitalizes Damaged Hair
- No Harmful Side Effects
- Drug-Free Treatment

**CLINICALLY PROVEN**

**FDA CLEARED\***  
FOR MEN & WOMEN

**DOCTOR RECOMMENDED**

## CLINICALLY PROVEN

HairMax has conducted more studies on light energy for hair growth than any other company to date. HairMax laser devices have been proven effective in 7 clinical studies at top medical research centers in the U.S.

Clinical study results have been published in 6 peer reviewed medical journals. HairMax laser devices have been granted 7 FDA Clearances and 14 medical device licenses world-wide.

over **90%** of clinical study participants experienced significant hair growth.+

**129** additional hairs per sq. inch, average increase in hair count.+

+ Based on clinical studies

\*HairMax laser devices are indicated to treat Androgenetic Alopecia, and promote hair growth in males who have Norwood Hamilton Classifications of IIa to V and in females who have Ludwig (Savin) I-4, II-1, II-2, or frontal patterns of hair loss and who both have Fitzpatrick Skin Types I to IV. Laser Light - Avoid direct eye exposure.

## HAIRMAX® LASERBAND 82

HairMax is the pioneer of at-home laser hair regrowth technology and the global leader in the field, bringing doctor recommended laser treatments to hair loss sufferers around the world.

The LaserBand 82 is a game changer in the hair loss industry, featuring 82 medical grade lasers, an innovative lightweight band design and the fastest treatment available – as little as 90 seconds, 3 days per week.

## SEEING IS BELIEVING




**Female:** Significant new hair growth in front temple area. Individual results may vary.



**Male:** Significant new hair growth, density and fullness. Individual results may vary.

## 3 ELEMENTS FOR EFFICACY

Our proprietary design includes 3 elements for efficacy – patented hair parting teeth, laser light intensity and laser light density – providing a highly effective hair loss treatment.



Hair is designed by nature as a photo protectant, its purpose is to block sunlight from reaching the scalp. The teeth on HairMax laser devices part your hair to help deliver the light energy to the follicles for optimal results.

